

Tests for molecular diagnostics of HLA associated diseases

IVD

CE

RNDr. David Elsnic

Tests for molecular diagnostics of HLA associated diseases



BAG Health Care GmbH develops and distribute worldwide wide range of HLA tests. Recently we focus on molecular tests and produce both SSP and SSO tests. Due our big experiences in HLA field we produce high quality tests for HLA associated diseases. We keep out tests up to date with latest knowledge about associations to provide for our customers most accurate diagnostics.

- HISTO TYPE Celiac Disease
- HISTO TYPE Narcolepsy
- HISTO TYPE B57
- HISTO TYPE B27

BAG HISTO SPOT DQB1 / DQA1

Celiac disease

Celiac disease is an autoimmune reaction triggered by gluten which is an ingredient of different cereals. If not diagnosed early this leads to chronic inflammation and destruction of the small intestine. Celiac disease is strongly associated with the DQA1*05:01- DQB1*02:01 and DQA1*03-DQB1*03:02 haplotype. Additionally, DR3, DR7 and DR11 alleles can be used as genetic markers.

BAG HISTO TYPE COELIAC DISEASE HAPPY PACK





- Simple SSP method with products identification on gel.
- Interpretation either manually with diagram or with our dedicated software Histo Match..
- 23 PCR mixes designed for precise typing of all relevant alleles (HLA-DQ2 a HLA-DQ8) which are associated with CD.
- It is not full HLA typization, so therefore: HLA specialist not needed for interpretation.
- Ideal in combination with serological EIA tests.
 - HLA negative result is important (exclusion of CD)
 - Positive serological result is important (confirmation of CD)
- Non invasive method
- Fast and for reasonable price
- Results valid for whole life (EIA and biopsy has to be repeated).
- DRB1 region included for better results interpretation.
- Standard method in whole Europe.
- **Patients suspicious for CD.**
- **Relatives of CD patients.**
- **Confirmed CD patients for severances prediction** (2xDQ2>DQ8+DQ2 > DQ2 > DQ8)

BAG HISTO SPOT DQB1 / DQA1



- SSO PCR test for coeliac disease diagnose for users with high throughput of samples. With this test you can easily interpret up to 200 samples per day. Assay is performed automatically on MR. SPOT processor and interpreted also automatically with BAG software – HISTO MATCH.



Auswertediagramm / Evaluation diagram

HISTO TYPE Celiac Disease

<p>HLA DRB1 – DQA1 – DQB1 associations:</p> <p>1. DRB1*03 - DQA1*05:01 - DQB1*02:01 2. DRB1*04 - DQA1*03:01 - DQB1*03:02 3. DRB1*07 - DQA1*02:01 - DQB1*02:02 DRB1*11 - DQA1*05:05 - DQB1*03:01 4. DRB1*07 - DQA1*02:01 - DQB1*02:02</p>	<p>Reaction pattern:</p> <p>1, 5, 6, 19 2, 12, 13, 14, 15, 18 3, 4, 6, 7, 9, 10, 11, (15), 17, 23 3, 6, 7, 17</p>	<p>Susceptibility to Celiac Disease:</p> <p>high → high → very high } high low </p>
---	---	--

In case of a HLA-DQB1*02 homozygous sample, the risk of CD development is 5 times higher than by a DQB1*02 heterozygous sample

Positive Reactions		Lenght of Internal control in bp																									
Spezifitäten / Specificities		1070	1070	1070	1070	429	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070
Ser.Type		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		
DR3	DRB1*03:01:01:01:01:03,01:05-41,43-67N,68N-75,77-80,82,83 / *11:07,53,103,105,107,125 / *14:104,111w / DRB3*01:14	1																									
DR4	DRB1*04:01:01:05:11w,05:12-62,64-81N,82-94N,95-109,111-119N,120N-124 / *03:49 / *08:38 / *11:102:01,102:02		2																								
DR7	DRB1*07:01:01:01:01:01,01:01:02,01:02w-01:04,03-09,10N-21			3																							
DR11	DRB1*11:01:01:01:01:21,23-29,31w,32,33w,34,35w-37:02,39-52w,53-65:02w,66,68-70,72-96w,97-129 / *03:08,65				4																						
DQ2	DQB1*02:01:01:01:01:05,04,07?					5!	6																				
DQ2	DQB1*02:02,06w						6	7																			
DQ2	DQB1*02:03							7	8																		
-	DQB1*02:05					?	6	?																			
DQ7(3) / -	DQB1*03:01:01:01:01:01:03,01:03-01:06,16,19,21,22,24,27-29,35,36,42									9	10	11					15										
DQ7(3)	DQB1*03:01:02									9	10	11															
DQ8(3) / -	DQB1*03:02:01,02:02,02:04,05:03,05:04,11,18,32,37													12	13	14	15										
DQ8(3)/DQ3	DQB1*03:02:03,05:01,05:02													12		14											
DQ9(3) / -	DQB1*03:03:02:01-03:02:03,03:04,15,17,26,30,31,33,34,38,39,41,43																14	15									
DQ7(3)	DQB1*03:04									9	10	11	12	13		15	16										
DQ3 / DQ4	DQB1*03:03:03,06,20,25 / DQB1*04:01:01-07																14										
-	DQB1*03:02:05,07													12	13	15											
-	DQB1*03:08													12	13	14											
-	DQB1*03:09,44									9	10						15										
DQ8(3)	DQB1*03:10										10	11					14	15									
-	DQB1*03:12										10						14	15									
-	DQB1*03:13																15										
-	DQB1*03:14										10		12	13		15	16										
DQ4	DQB1*04:02:02								8								14										
-	DQB1*04:08											11				14											
-	DQA1*02:01																	17									
-	DQA1*03:01:01																			18							
-	DQA1*05:01:01:01,01:01:02,01:02?																				19						
-	DQA1*05:02																					20					
-	DQA1*05:03,06,07																				19		21*				
-	DQA1*05:04																				?			22			
-	DQA1*05:05:01:01-01:01:03,08-11																									23	
-	Kontaminationskontrolle/Contamination Control																										24
	Länge in bp / length in bp	220	200	140	177/171	800!	150	140	150	105	200	185	120	95	145	135	135	105	485	235	90	200	205	235			

w = schwache Reaktion / weak reaction ! / ! Die Größe der internen Kontrolle in Mix 5 ist 429bp / The size of the internal control in mix 5 is 429bp !

* Die spezifische Bande in Mix 21 kann schwächer sein als die anderen spezifischen Reaktionen. / The specific band in Mix 21 may be weaker than the other specific bands.

Interpretation with Histo Match software




- Dedicated software HISTO MATCH for our SSP and SSO test if for our customers free of charge.
- Suitable for:
 - easy and fast interpretations
 - Creating and printing reports
 - Protected results database
- HM for BAG HistoType Celiac disease kit:
 - Easy interpretation with virtual gel picture
 - Missing controls tests interpretation
 - Clean typization reports
 - Precise model for CD predisposition
 - Possible rare alleles included in report
 - Results database

Patient:
Sample ID 720417 Status Initial

Celiac Disease DQB1_316S1 (IMGT v. 3.13.1)

1

All 4D CWD User



Coeliac

Summary **DQ2.x/DQ2 [Coeliac Risk: High]**

DRB1*07:01:01:01/07:01:01:02/07:01:03/07:01:04/07:03/07:04/07:05/07:06/07:07/07:08/07:09/07:10N/07:11/07:12/07:13/07:14
 DQB1*02:02:01/02:02:02/02:10/02:11
 DQA1*02:01

HLA-DRB1
Suggested Ambiguous Result
DRB1*07,-; DRB1*07,
DRB1*07


HLA-DQB1
Suggested Ambiguous Result
DQB1*02,-; DQB1*02,
DQB1*02; DQB1*02,
DQB1*03

HLA-DQA1
Suggested Ambiguous Result
DQA1*02,-; DQA1*02,
DQA1*05

Results History Comments

Hit Table DQB1_316S1 IMGT 3.13.1 Version 18.7.2013

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
	220	200	140	175	800	150	140	150	105	200	185	120	95	145	135	135	105	485	235	90	200	205	235
Group/Allele	MM	CWD																					
DQB1*02,-(2)	0	[2]																					
DRB1*07, DRB1*	0	[1]																					
DQB1*02, DQB1*	0	[1]																					



Patient: **Celiac Disease** DQB1_316S1 (IMGT v. 3.13.1)

Sample ID 720417 Status Initial

Buttons: [Print] [1] [Cancel] [OK] [Checkmark]

Buttons: All 4D CWD User

Coeliac

Summary

Coeliac Risk: High

The diagram shows two HLA-DQ alleles. On the left, the DQ2.x allele is shown with a purple β-chain (DQB1*02) and an orange α-chain (DQA1*x). On the right, the DQ2 allele is shown with a green β-chain (DQB1*02) and a white α-chain (DQA1*x).

- HLA-DRB1**
Suggested Ambiguous Result
DRB1*07,-; DRB1*07,
DRB1*07
- HLA-DQB1**
Suggested Ambiguous Result
DQB1*02,-; DQB1*02,
DQB1*02; DQB1*02,
DQB1*03
- HLA-DQA1**
Suggested Ambiguous Result
DQA1*02,-; DQA1*02,
DQA1*05

Buttons: Results History Comments

Buttons: [Checkmark]

Hit Table DQB1_316S1 IMGT 3.13.1 Version 18.7.2013

			1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
			220	200	140	175	800	150	140	150	105	200	185	120	95	145	135	135	105	485	235	90	200	205	235
Group/Allele	MM	CWD																							
▶ DQB1*02,-(2)	0	[2]																							
▶ DRB1*07, DRB1*	0	[1]																							
▶ DQB1*02, DQB1*	0	[1]																							

HISTO MATCH Patient Report



Patient: 720417 28.1.2014

Name: David elsnic

Gender: Male

Blood Group: A + CMV: -

Comments:

Sample ID: 720417 Collected: 28.1.2014

Sample Received: 28.1.2014 16:33:37

Extraction Method:

Tests

Tests

Worklist

Extraction

Date Performed

Analyzed By david

Date Analyzed 28.1.2014

Approved By

Date Approved

Authorized By

Date Authorized

Kit

Batch: DQB1_31651

Expiry: 30.9.2015

Hit Table: DQB1_31651, IMGT db: 3.13.1, ver:18.7.2013

Method: SSP

Suggested Ambiguous Result

Summary: DQ2.x/DQ2 [Coeliac Risk: High]

Filter: Display All Alleles

Serology: DR7/DRNull/DR-;DQ2/DQ-;DQA1?

Assigned NMDP: DRB1*07:JKZW, DRB1*07:JKZW;DQB1*02:BPPD, DQB1*XX:XX;DQA1*02:01, DQA1*XX:XX

Possible NMDP: DRB1*07:JKZW, DRB1*07:JKZW;DQB1*02:BPPD, DQB1*XX:XX;DQA1*02:01, DQA1*XX:XX

Group Result:

DRB1*07,- [2]

DRB1*07, DRB1*07 [1]

DQB1*02,- [2]

DQB1*02, DQB1*02 [1]

DQB1*02, DQB1*03 [1]

DQA1*02,- [2]

HISTO MATCH Patient Report



Patient: 720417 28.1.2014

DQA1*02, DQA1*05 [1]

Allelic Result:

DRB1*07:01:01:01/ 07:01:01:02/ 07:01:03/ 07:01:04/ 07:03/ 07:04/ 07:05/ 07:06/ 07:07/ 07:08/ 07:09/ 07:10N/ 07:11/ 07:12/ 07:13/ 07:14/ 07:15/ 07:16/ 07:17/ 07:18/ 07:19/ 07:20/ 07:21, DRB1*07:01:01:01/ 07:01:01:02/ 07:01:03/ 07:01:04/ 07:03/ 07:04/ 07:05/ 07:06/ 07:07/ 07:08/ 07:09/ 07:10N/ 07:11/ 07:12/ 07:13/ 07:14/ 07:15/ 07:16/ 07:17/ 07:18/ 07:19/ 07:20/ 07:21

DRB1*07:01:01:01/ 07:01:01:02/ 07:01:03/ 07:01:04/ 07:03/ 07:04/ 07:05/ 07:06/ 07:07/ 07:08/ 07:09/ 07:10N/ 07:11/ 07:12/ 07:13/ 07:14/ 07:15/ 07:16/ 07:17/ 07:18/ 07:19/ 07:20/ 07:21, DRB1*07:01:02

DQB1*02:05, DQB1*02:05

DQB1*02:02:01/ 02:02:02/ 02:10/ 02:11, DQB1*02:02:01/ 02:02:02/ 02:10/ 02:11

DQB1*02:05, DQB1*02:07/ 02:08/ 02:09

DQB1*02:05, DQB1*02:02:01/ 02:02:02/ 02:10/ 02:11

DQB1*02:05, DQB1*02:06

DQB1*02:07/ 02:08/ 02:09, DQB1*02:02:01/ 02:02:02/ 02:10/ 02:11

DQB1*02:02:01/ 02:02:02/ 02:10/ 02:11, DQB1*02:06

DQB1*02:05, DQB1*03:23

DQB1*02:02:01/ 02:02:02/ 02:10/ 02:11, DQB1*03:23

DQA1*02:01, DQA1*02:01

DQA1*02:01, DQA1*05:01:02

Could Not Exclude

Positive Reactions

3, 6, 7, 17

Ignored Reactions

Comments

NARCOLEPSY

- Narcolepsy belongs to a groups of sleeping disorders (Dyssomnias) with symptoms like daytime somnolence, sleep paralysis and hallucinations. 98 % of the Caucasian patients with narcolepsy have the DQB1*061:02 haplotype. Therefore, HLA typing is helpful to confirm or exclude a diagnosis.
- Prevalence is 0,02–0,18% of whole population, (depending on ethnicity).
- Association between narcolepsy and HLA is long known.
- 85-95% narcolepsy patients carry HLA haplotype DRB*15:01-DQA1*01:02-**DQB1*06:02**, in common population is this haplotype present just in 30%
- Strongest association is towards subtype **DQB1*06:02**, which therefore constitutes best marker for narcolepsy.
- *Test great for exclusion (when results negative) of narcolepsy, positive results has weak prediction value.*

BAG HISTO TYPE Narcolepsy



- Simple SSP method with amplicons interpretation on gel
- Interpretation manually with diagram
- Just 8 PCR mixes for detailed examination of all HLA-DR, DQA and DQB alleles known to be associated with narcolepsy.
- It is not full HLA typization, just selected alleles, therefore HLA specialist not required for interpretation
- Non invasive method
- Fast and for reasonable price
- Results valid for whole life
- *All patients suspected for narcolepsy.*

Abacavir

- Treatment with antiretroviral drugs containing the active substance Abacavir (mainly in HIV therapy) is only allowed when the existence of HLA-B*5701 is excluded by molecular genetic methods, since otherwise it is possible that a hypersensitivity reaction against Abacavir may occur.

HISTO TYPE B57 Happy Pack

- SSP PCR kit for exclusion of HLA-B*57:01. Contains 7 mixes + contamination control.

Morbus Bechterew et al.

- There is a strong association between HLA-B27 and a group of seronegative arthritic diseases like Morbus Bechterew or Morbus Reiter. Determination of HLA-B27 is an important tool in the diagnosis of these diseases.

HISTO TYPE B27 Happy Pack

- SSP PCR kit for detection of HLA-B*27. Contains 1 mix.

HISTO TRAY B27

- Srological test for B*27 detections.